Mystery Powder

Name

Mystery Chemicals





Do not taste any of the powders. Goggles must be worn during this investigation. Iodine will cause a permanent stain. Tie back hair and loose clothing around an open flame.

Background Information: Chemists are scientists who study what different chemicals are shaped like, what color they are, how they behave, and how they react with other chemicals. What color or shape a chemical has is called a **physical property**. How a chemical behaves and how it reacts with other chemicals is called a **chemical property**. Chemists use both chemical and physical properties to identify and separate different kinds of chemicals.

For some chemists and forensic scientists, the identification and analysis of unknown substances is a daily task. For example, when a law enforcement officer discovers a suspiciously concealed white powder at a crime scene, it is the chemist's job to figure out whether the substance is cocaine, heroin, or just plain old table sugar. Using scientific problem-solving skills, the unknown can be correctly identified in most cases. This is known as **qualitative analysis**.

Problem: To identify 5 unknown samples of matter by their physical and chemical properties (qualitative analysis).

5 unknown white chemicals	lodine	Acid droppers
Spot Plate	Beaker with Water	Test tubes
Scoopulas	Eye droppers	Test tube rack
Hot Plate	Foil	

Materials

Procedure:

- 1. Analyze the chemical in the petri dish and record the appearance in column 2 in the observation chart.
- 2. Measure out 10mL of water in a small graduated cylinder.

Add this water to a test tube. Using a scoopula, add a small amount of unknown powder. Stopper and shake the test tube. Observe the solubility and record in column 3.

- 3. Using a scoopula, add a small amount of unknown powder to two separated wells in a spot plate.
- 4. Use an eyedropper to put 3 drops of acid in the first well. Observe and record observations in column 4.
- 5. Use a different eyedropper and put 3 drops of iodine solution (stains the skin) in the second well. Observe and record in column 5.
- 6. Use the foil to make a small boat. Put a small amount of unknown powder in the foil. Turn the hot plate on high (8-10). Put it on the hot plate and wait 1-2 minutes to see if a reaction occurs. Observe and record in column 6.
- 7. Rinse your spot plate completely and repeat steps 1-6 with each unknown powder. Wash and dry the spot plate between each use. Use a new piece of foil for each powder.

Da	ata:				
#1	#2	#3	#4	#5	#6
	Appearance	Solubility in	Acid	lodine	Heat
POWDER	(colour, crystal	Water (soluble	(bubbles/fizz or	(brown or black)	(describe the
	form, state)	or insoluble)	no reaction)	(reaction or
	, ,		norodotiony		no reaction)
					noreautiony
1					
2					
3					
Ŭ					
Λ					
-					
E					
5					

Data Analysis: Use the Properties Chart to compare your observations. Identify each unknown powder:

1.	
2.	
3.	
4.	
5.	

Questions: (5)

1. What properties did **all of the chemicals** have in common?

2. How did you tell the difference between powders 1 and 2 (list minimum 2 properties to distinguish between them)?

3. How did you tell the difference between powders 3 and 4 (list minimum 2 properties to distinguish between them)?

Thinking / Inquiry: _____/20